



Versatile Scientific Products
A Division of.....
Grizzly Analytical



PO Box 573

E. Setauket, New York, 11733

Main: (631) 513-3300

Tele: (631) 444-5616

FAX: (631) 444-5627

info@VersatileScientificProducts.com

www.VersatileScientificProducts.com

TECHNICAL DATA SHEET

SHOE COVERS

Product Code: APP0330-LTX

Colors Available: White

Sizes Available: Universal Extra Large

Packing: Poly Bag
100 pieces per bag
3 bags/ 300 pieces per case

Characteristics

White Compressed Polyethylene (CPE), commonly called Polylatex, is made from a mixture of 65% Low Density Polyethylene (LDPE) and 35% Light Low Density Polyethylene (LLDPE). It is practically impervious and is excellent for containing moisture and dirt. Polylatex is non-toxic and odorless. Polylatex has excellent chemical resistance and great strength. The material is heavily textured for better grip.

Specifications

Material:	Polylatex	
Dimension:	<u>Universal</u>	<u>Extra Large</u>
a) Length:	13"	18"
b) Width:	6 ½"	8"

Product Details

- Elastic at the ankle provides a secure, comfortable closure.
- Poly Latx film has anti-static properties
- Lint-Free Material