



Versatile Scientific Products
A Division of.....

Grizzly Analytical 

PO Box 573

E. Setauket, New York, 11733

Main: (631) 513-3300

Tele: (631) 444-5616

FAX: (631) 444-5627

info@VersatileScientificProducts.com

www.VersatileScientificProducts.com

MATERIAL SAFETY DATA SHEET

POWDER FREE NITRILE EXAM GLOVES

Identification

Product Code: GLP0065

Product Description: Non-Sterile, Ambidextrous, Disposable, Powder-Free Nitrile Gloves, Blue Colored.

CASE #: Not Established

IATA/IACO Hazard Rating: Not Restricted

D.O.T Rating: Not Restricted

Physical Data

Composition: Nitrile Potassium Hydroxide
Sulfur Titanium Dioxide
Accelerator, ZDBC Stabilizer
Zinc Oxide

Fire and Explosion

Stability: Stable

Flash Point: N/A

Extinguishing Media: Foam
Carbon Dioxide
Water

Incompatible Materials: None

Normal Precautions

Ventilation: Not necessary under normal conditions

Equipment: Not necessary

Spills Clean-up: N/A

Exposure Data

Eyes: N/A

Storage: Keep away from sunlight and heat. Store at room temperature, 10-40°C.