



**Versatile Scientific Products**  
A Division of.....

**Grizzly Analytical**



**PO Box 573**

**E. Setauket, New York, 11733**

**Main: (631) 513-3300**

**Tele: (631) 444-5616**

**FAX: (631) 444-5627**

**info@VersatileScientificProducts.com**

**www.VersatileScientificProducts.com**

## **MATERIAL SAFETY DATA SHEET**

### **POWDER FREE VINYL GLOVES**

#### **Identification**

Product Code: GLP0080      GLP0081      GLP0085

Product Description: Powder Free Vinyl Glove

CASE #: Not Established

IATA/IACO Hazard Rating: Not Restricted

#### **Physical Data**

Composition: Synthetic PVC - Polyvinyl Chloride Polymer  
Plasticizer DINP  
Plasticizer Epoxidized Soybean Oil  
Plasticizer TXIB  
Stabilizer Calcium-Zinc

#### **Fire and Explosion**

Stability: Stable

Flash Point: N/A

Extinguishing Media: Foam  
Carbon Dioxide  
Water

Conditions to Avoid: None

Incompatible Materials: None

#### **Normal Precautions**

Ventilation: Not necessary under normal conditions

Equipment: Not necessary

Spills Clean-up: N/A

#### **Exposure Data**

Eyes: N/A

#### **Other Precautions**

Storage: No special precautions necessary. Keep at room temperature.