Nitrile Examination Gloves







Confidence is key. When you touch, smell and stretch a pair of our Nitrile Gloves, you will feel the quality and consistency that is inherent in all **DuraShield** gloves. Slip your hand into a pair, and you will quickly discover that these are gloves you can rely on for protection and comfort.

Precision-Grip[™] **Nitrile Gloves** are entirely LATEX-FREE making them suitable for latex sensitive users. Strong, yet soft and flexible, these gloves are comfortable for extended wear purposes. Precision-Grip Nitrile offer unsurpassed tactile sensitivity, dexterity, flexibility and tensile strength.

Available in Lightly Powdered or Powder Free, with a beaded cuff for easy donning. Ambidextrous gloves are nine inches in length with micro-textured fingers for secure wet/dry handling. **Precision-Grip** Nitrile Exam Gloves are the ideal choice for medical and laboratory use as well as industries where hand protection is required.

Specifications

Tensile Properties

<u>Puncture Resistance</u>	Not Le	ss than	5.0 N
<u>Unaged</u>			
Tensile Strength (MPA)	Not Le	ss than	14.0
Elongation, %	Not Les	ss than	500
Force at Break	Not Le	ss than	7.5 N
Aged 70 +/- 2° C for 7 days			
Tensile Strength (MPA)	Not Le	ss than	14.0
Elongation, %	Not Les	ss than	500
Force at Break	Not Le	ss than	5.5 N

*Specifications subject to change without notice

Dimensions

<u>Width</u>
Small
Medium
Large 106 ½ 5 mm
Extra Large 116 ½5 mm
<u>Length</u> Not Less than 240 mm
Thickness (single wall) Not Less than 0.08 mm



Nitrile Examination Gloves





Precision-Grip LP Lightly Powdered Nitrile Gloves

	3 3			
Product No.	Description	Size	QTY/BOX	QTY/CS
GLP0055-S	Lightly Powdered Nitrile Exam Gloves – 9"	Small	100	1,000
GLP0055-M	Lightly Powdered Nitrile Exam Gloves – 9"	Medium	100	1,000
GLP0055-L	Lightly Powdered Nitrile Exam Gloves – 9"	Large	100	1,000
GLP0055-XL	Lightly Powdered Nitrile Exam Gloves – 9"	X-Large	100	1,000

Precision-Grip PF Powder Free Nitrile Gloves

Product No.	Description	Size	QTY/BOX	QTY/CS
GLP0065-S	Powder Free Nitrile Exam Gloves – 9"	Small	100	1,000
GLP0065-M	Powder Free Nitrile Exam Gloves- 9"	Medium	100	1,000
GLP0065-L	Powder Free Nitrile Exam Gloves- 9"	Large	100	1,000
GLP0065-XL	Powder Free Nitrile Exam Gloves- 9"	X-Large	100	1,000



Additional

Gloves Include:

Vinyl Exam Powder-Free GLP0081 (S-XL)
Vinyl Exam 12" Powder Free GLP0088 (S-XL)
Vinyl Exam Lightly Powdered GLP0091 (S-XL)
Latex Exam Powder-Free COAP110 (XS-XL)
Latex Exam Lightly Powdered COAP110 (XS-XL)
Latex Exam Powder-Free HIGH RISK 12" GLP0110 (S-XL)
CLEANROOM GLOVES ARE ALSO AVAILABLE

7intri	hutad	D. /.
ハシロ	buted	BA:

