

# Borosilicate Glass Culture Tubes

Globe Scientific offers a high quality line of borosilicate glass culture tubes. These tubes offer exceptional quality at great prices. They are neatly packaged in heavy duty boxes of 250 and are shrink wrapped to keep them particle free. The boxes are compact and designed to fit into most test tube dispensers.



- Made from superior quality borosilicate glass to reduce pH changes and contaminants that are potentially leached from soda lime glass tubes
- Tubes feature well-rounded bottoms and smooth, fire-polished rims
- Dimensional specifications are rigidly controlled for reliable, uniform results
- Produced in an ISO 9001:2008 certified factory
- Tubes can be used for tissue culture and blood banks



**TUBESAVER™  
PACKAGING SYSTEM!**  
Our **heavy duty** corrugated boxes prevent **breakage** during shipping!



Item #	Size	Overflow Capacity	Packaging	Unit
1503	10 x 75mm	4mL	250/box, 4 boxes/case	1000
1505	12 x 75mm	6mL	250/box, 4 boxes/case	1000
1510	13 x 100mm	10mL	250/box, 4 boxes/case	1000
1512	16 x 100mm	14mL	250/box, 4 boxes/case	1000
1515	16 x 125mm	19mL	250/box, 4 boxes/case	1000
1517	16 x 150mm	23mL	250/box, 4 boxes/case	1000
1519	18 x 150mm	29mL	125/box, 4 boxes/case	500
1522	20 x 150mm	36mL	125/box, 4 boxes/case	500



# Diamond™ Culture Tube Caps

These popular caps fit plain end culture tubes and flasks with DeLong style necks or ports. The caps are molded from durable, autoclavable polypropylene (PP). They are available in a full range of five sizes and five colors for easy sample identification and organization. Recommended tube diameter is indicated on each cap.

- Polypropylene (PP)
- Autoclavable
- Ridges inside cap provide secure grip and gas exchange
- Available in five sizes and five colors
- Maintain sterility without cotton plugging
- Reduce sample evaporation compared to cotton plugging
- Provide same evaporation and sterility benefits as stainless steel closures

**By popular demand!**



Max Tube Diameter	Bag Qty.	Unit	Item#				
			Natural	Blue	Green	Red	Yellow
13mm	100	1000	118150C	118150B	118150G	118150R	118150Y
16mm	100	1000	118152C	118152B	118152G	118152R	118152Y
18mm	100	1000	118154C	118154B	118154G	118154R	118154Y
20mm	100	500	118156C	118156B	118156G	118156R	118156Y
25mm	100	500	118158C	118158B	118158G	118158R	118158Y

## Dispensers for Glass Culture Tubes

These convenient dispensers are made from heavy duty white painted metal. Each holds 1 box of 250 tubes. The front acrylic cover pulls up to dispense the tubes.

- Conserve valuable benchtop real estate
- Keep tubes clean until point of use
- Safely remove only what you need

Item #	Description	Unit
1590	For 10 x 75mm and 12 x 75mm tubes	Each
1592	For 13 x 100mm tubes	Each



[www.globescientific.com](http://www.globescientific.com)

**Customer service, sales and technical support:**

Phone 1 (800) 394-4562 • 1 (201) 599-1400

Fax 1 (201) 599-1406

E-mail [mail@globescientific.com](mailto:mail@globescientific.com)